

Quantitative Analysis of C-Phycocyanin from *Spirulina Pacifica* (Low Temperature Method)

Background:

Spirulina and other blue-green algae contain c-phycoyanin, which acts as an accessory pigment when light energy is captured and transferred to chlorophyll a. This is a spectrophotometry method adapted to extract and quantify a relatively pure c-phycoyanin fraction from *Spirulina Pacifica*.

Equipment and instruments:

Spectrophotometer at 620nm
Refrigerator (4 C)
Phosphate buffer (pH 7.0)*
10 ml centrifuge tubes
Cooled centrifuge (10 C @ 3500 RPM)
Dessicator
Weigh pans
Analytical balance

Methods:

Dry Weight

- 1) Place drying pans in oven for 30 minutes place in desiccator to remove excess moisture.
- 2) When pans are cool, weigh and record weight of pan.
- 3) Tare the balance with the pan on it and place about two grams of powder in the pan.
- 4) Record the weight of the powder.
- 5) Place pan and powder in the oven to dry for six hours.
- 6) Remove pan and powder from the oven and place in desiccator 15 minutes to cool.
- 7) Weigh and record the total weight of the pan and the dry powder.
- 8) Perform duplicates for each sample.

Dry Weight Calculations

$$\text{Percent dry wt} = \frac{\text{pan (g)} + \text{dried powder (g)} - \text{pan wt (g)}}{\text{powder wt (not dried) (g)}}$$

Phycocyanin Assay

- 1) Weigh accurately 40 mg. *Spirulina* powder into a 10-ml. centrifuge tube.

- 2) Add 10 mls. of the 100 mM phosphate buffer (100-mM Phosphate buffer contains 10.64 g. K₂HPO₄ and 5.29g. KH₂PO₄ per liter, pH 7.).
- 3) Vortex to mix well.
- 4) Store in refrigerator overnight.
- 5) Vortex to mix well.
- 6) Centrifuge 5 minutes at 10 C at 3500 RPM.
- 7) Read absorbency of each replicate at 620 nm, using phosphate buffer as blank.
- 8) Average absorbency readings for dilution replicates.
- 9) Calculate percent C-Phycocyanin:

Derivation of pure C-Phycocyanin:

$$\% \text{ pure CPC} = \frac{A_{620} \times (10) \times (100)}{7.3 \times (\text{mg. sample}) \times (\% \text{ dry wt.})}$$

where 7.3 is Extinction coefficient of CPC at 620 nm
 10 is total volume;
 100 represents 100%.

Derivation of crude C-Phycocyanin:

$$\% \text{ crude CPC} = \frac{A_{620} \times (10) \times (100)}{3.39 \times (\text{mg. sample}) \times (\% \text{ dry wt.})}$$

where 3.39 is Extinction coefficient of CPC at 620 nm
 10 is total volume;
 100 represents 100%.

References

Boussiba S. and A. Richmond. 1979. Isolation and purification of phycocyanins from the blue-green alga *Spirulina platensis*. Arch. Microbiol. 120:155-159.

Spirulina Pacifica Technical Bulletin #004
 Contact: Dr. R. Todd Lorenz
 Cyanotech Corporation
 Phone: 808-326-1353
 FAX: 808-329-3597
 Email: florenz@kona.net
www.cyanotech.com

Revision Date: November 24, 1998